Instruction Manual

OMRON

Body Composition Monitor


Model **HBF-375**Karada Scan

■ Illustration in the book is for reference only.

Contents

How to use HBF-375	1
Foreword Notes on SafetyKnow Your Unit	
Before Measuring (various settings Installing the Battery Battery Life and Replacement Setting Date and Time. Setting Personal Data (Date of Birth, Gender, Height) Changing Personal Data Deleting Personal Data	. 8 8 9 11
Measuring Methods and Measuring Results Measuring Weight and Body Composition (with Personal Profile Number Button / GUEST Button) Using the Memory Function	15 22 24 25
Care and Maintenance Error Displays	27 27 28
Appendix – Application Guide Information on Body Composition Same Age Comparison Feature and Measuring Principle of Body Composition Monitor Things to Take Note Q&A of body composition monitor	35 36 37
_	

- Before using this unit for the first time, please be sure to read this Instruction Manual carefully and use the unit safely and properly.
- Please keep this Instruction Manual at hand all the time for future reference.

How to use HBF-375

Before measuring (various settings)

Please make the following preparation before using for correct measuring.

Installing the battery (page 8)

2 Setting date and time (P page 9)

3 Setting personal data

Measuring methods and measuring results (© page 15)

1 Connecting the power

Turn on the power.

Do not take out display unit before "0.0 kg" is displayed.

2 Selecting personal profile number

Press Personal Profile Number button to select "personal profile number".

Please set personal data before initial use. (Fee page 11)

In addition to setting personal data, GUEST mode may also be used to measure weight, body fat and skeletal muscle, etc. (P page 16)

3 Stepping on the main unit

Take a measurement bare-footed.

Use correct posture for measurement. (Fig. page 17 and 18)


Measuring...


Measuring completed

Measure on flat and hard ground

Notes on Safety

Please read carefully before use.

- The following is intended to prevent personal injury or property damage of the user or other people caused by the product.
- Symbols and definitions are as follows.

Meaning of danger, warning and caution

<u> </u>	Improper use may cause danger resulting in death or serious injury.
⚠ Warning	Improper use may result in possible death or serious injury.
⚠ Caution	Improper use may result in injury or property damage*.

* Property damage includes damage to house, property, livestock and pet.

Meanings

0	
Compulsory	

Compulsory item (which must be observed).

○ Forbidden item (which is not permitted).

⚠ Danger

When using

Never use this unit in combination with the following medical electronic devices.

- 1 Medical electronic implants such as pacemakers.
- 2 Electronic life support systems such as an artificial heart/lung.
- ③ Portable electronic medical devices such as an electro-cardiograph.
- This unit could cause these devices to malfunction, posing a considerable health risk to users of these devices.

A Warning

When using

Never start weight reduction or exercise therapy without the instructions of a doctor or specialist, and do not judge by yourself.

• If you judge by yourself, it might cause damage to your health.

Marning

Installation and usage

Do not use the unit on slippery surfaces, such as ceramic tile or wet floor.

• Otherwise, it may cause you to slip and result in injury.

Measurement

People with disabilities should always be assisted by another person when using this unit.

• Otherwise, you might slip. Take care when measuring.

Compulsory

Do not use this unit when your body and / or feet are wet, such as after taking a bath.

- Otherwise, it may cause you to slip and result in injury.
- Otherwise, water will come into the unit to cause fault. Correct measurement may be impossible.

Do not jump onto the unit, or hop up and down on the unit.

• Otherwise, it might cause slip and injury. Impact exceeding allowed range may cause damage.

Do not step on the edge or display area of the main unit.

 Otherwise, it might cause slip and injury and correct measurement may be impossible.

Forbidae

When using

This unit is not intended for professional use (for example, hospitals or other medical facilities).

 This unit is intended for home use only, it does not support the standards required for professional use.

Do not use a cellular phone near the display unit or the main unit.

· Otherwise, it might cause abnormality.

Application scope

The unit is designed to measure weight and resistance and obtain BMI, body fat and skeletal muscle, etc. through analysis of weight and resistance.

M BMI, also known as body weight index, is a physique indicate index. BMI = weight (kg) / square of height (m).

♠ Caution

Installation and usage

Do not place the unit against the wall.

• Otherwise, it might cause slip, injury or fault.

Forbidden

Remove the display unit from main unit before stepping on the main unit.

• If you try to remove the display unit while stepping on to the main unit, you might slip and cause injury due to lost balance.

When inserting display unit into the main unit, do not touch inside part of grip electrodes metal.

· Otherwise, it might hurt your finger.

Forbidden

Do not disassemble, repair, or remodel the display unit or the main unit.

· Otherwise, it might cause injury or fault.

Forbidden

Measurement

Taking a measurement barefooted.

• Otherwise, it may cause you to slip and result in injury. And it will affect correct measuring of body fat and other index.

Battery usages

Do not insert the batteries with the (+)(-) polarities in the wrong direction.

• Otherwise, it might leak, heat or break and cause you injury.

Compulsory

Do not use new and worn batteries and different types of batteries together.

• Otherwise, it might leak, heat or break and cause you injury.

Compulsory

Remove the batteries from this unit when you are not going to use it for a long period of time (approximately three months or more).

Use batteries specified for this unit.

• Otherwise it might leak, damage the main unit or cause person injury.

Compulsory

Storage

Keep the unit out of the reach of young children.

• Otherwise it might cause injury due to winding wire or falling down.

Notes

■ Notes on Installation

Keep the unit away from wet or watery place, direct sunlight, direct air conditioner wind and fire.

· Otherwise, it might cause fault.

Do not use the unit on grass mat, carpet or other soft ground.

 Otherwise, when stepping on the main unit, it will sunken and measurement will be affected.

Notes on use

Clean the unit before use if it is used by people whose hand or foot sicken for skin disease or other infective disease.

 Otherwise, it will cause infection. If use the unit together, use a cloth moistened with water or neutral detergent to squeeze it well before wiping the unit, then wipe dry with a dry cloth.

Do not use the unit for purpose other than measuring weight and body composition.

Do not draw the cord between the main unit and display unit with strength.

• Otherwise, the cord might come off and cause fault.

Since the product is precise instrument, please do not drop, violate or shock it.

· Otherwise, it might cause fault.

Do not move the unit by holding the display unit.

• Otherwise, the main unit might drop and be damaged or cause injury to personnel.

■ Notes on daily maintenance

Always keep the unit clean before use.

- Wipe the main unit with a soft dry cloth.
- If necessary, use a cloth moistened with water or neutral detergent and squeeze it well before wiping the unit, then wipe dry with a dry cloth.

Do not wash display unit or the main unit.

• Otherwise, it might cause fault.

Do not use volatile detergent or thinner, to clean the unit

 Otherwise it might cause color change or fault.

■ Notes on storage

Do not keep the product in the following places.


- · Watery place
- Place subject to high temperature, humidity, direct sunlight, dust or saliferous air.
- Places with the risk of inclination, vibration or sudden shock
- Places where chemicals are stored or where corrosive gas is present.

General Safety Precautions

- The unit might not correctly measure body composition for the following personnel.


 Children in growth age / people of high age / women of amenorrhea / patients having fever due to cold / pregnant women / patients having osteoporosis / patients having edema / patients receiving hemodialysis / professional sportsmen.
 - · Water content in their body and body composition might be very different from average value.
- For users under 10 years old, the unit displays measuring weight only.
- This unit cannot display visceral fat level, body age and same age comparison of users under 18 years old.
- This unit may measure body composition of users over 80 years old but only provide reference value.

Know Your Unit


Installing the Battery

Install four supplied AA batteries.

Turn over the main unit Open battery cover

- 1 Press the tab of the battery cover along arrow direction
- 2 Pull up the tab of the battery cover and open battery cover

Install the batteries in correct polarity as marked inside the battery compartment

The side with spring is \bigcirc .

3

Close the battery cover

■ Setting date and time (page 9)


Reset date and time after replacing the battery.

Battery Life and Replacement

- Four AA batteries will last approximately 1 year. (At room temperature 23°C, when measured 4 times per day.)
 Because the supplied batteries are for trial use only, they may have a shorter life.
- ★ When the low battery indicator flashes, it indicates running out of the battery.

Replace all four batteries with new ones (same type) at the same time.

- Replace the batteries after turning off the main unit power.
 - Personal data stored in memory are retained even if the batteries are removed.
 - Disposal of used batteries should be carried out in accordance with the city environment protection regulations for the disposal of batteries.

Setting Date and Time

When the power is turned on for the first time after changing the battery, it is necessary to set the date and time. The unit will calculate age according to current date, the registered personal data and date of birth. If not set correctly, measurement results would not be correctly displayed. The previous measurement results might be lost.

Turn on the power

Set "year"

Press (NEM () BUSS) button to select "year" and press DISPLAY button to set.

- Setting range is 2012 ~ 2030.
- Press continuously to rapidly shift up or down.

After setting "year", "month" flashes.

Set "month" and "date"

Press (MEM DELET) button to select "month" and press (DISPLAY) button to set.

After setting "month", date flashes.

button to select "date" and press (DISPLAY) button to set.

After setting "date", "hour" flashes.

4

Set "hour" and "minute"

- Press () button to select "hour" and press () button to set.
 - Time is 24-hour system.

 After setting "hour", "minute" flashes.
- Press button to select "minute" and press star button to set.

After the settings for year, month, day, hour and minute are displayed in that sequence, the power is automatically turned off.

Setting is now completed.

- When error occurs during setting...
 - Turn on the power. If "year" flashed, please start from step 2 of "Setting date and time" (IF page 9).
- The power is automatically turned off during setting...

 Turn on the power to reset. If no operation for 3 minutes during setting, the power will be automatically turned off.
- If you want to modify date and time which have already been set...

 Turn on the power. After "0.0 kg" is displayed, hold DISPLAY button more than 2 seconds. If "year" flashes, start from step 2 of "Setting date and time".

Setting Personal Data

(Date of Birth, Gender, Height)

For the measurement of body composition, it is necessary to set your personal data (age, gender, height). The unit is designed to store measurement results in the memory for 4 people.

* If no input of personal data in advance, you can use "GUEST" button to take measurement. (🕼 page 16)

1 Sentition Sent

With the display unit inside the main unit,

Turn on the power

After the display of current "year", "month and date" (MM/DD), it will display "0.0 kg". Given example of date "725" means July 25th.

Note: Do not step on the main unit or move the main unit before "0.0kg" is displayed. Do not put something on the main unit or touch display unit. Otherwise, it will display "Err (error)". (

| Page 27)

2 Omnon 3

Press the personal profile number button to select a personal profile number.

After setting personal profile number, "year of birth" will flash.

Note: When "0.0 kg" is displayed, it means that personal profile number has been registered successfully. If you want to use this number, please delete this personal data.

"Deleting personal data" (propage 14)

3

Set "date of birth"

- Press button to select "year" and press button to set.
 - Setting range is 1900 ~ 2030.
 - Press continuously to rapidly shift up or down.


After setting "year of birth", "month" flashes.

Press button to select "month of birth" and press obsplay button to set.

After setting "month of birth", "date" flashes.

After setting "date of birth", "male / female" flashes.

4

Set "gender"

Press button to select "male" / "female" and press (DISPLAY) button to set.

After setting "gender", "height" flashes.

5

Set "height"

Press (MEM) DUST) button to select "height" and press (DISPLAY) button to set.

• Setting range is 100.0 cm ~ 199.5 cm.

After the settings for age, gender, and height are displayed in that sequence, then 0.0 kg appears on the display.

Setting is completed.

You can now take a measurement.

■ When starting measurement...

Taking a Measurement (P page 17)

■ When setting personal data of other people...

Setting Personal Data (page 11)

■ After operation is completed...

Press power switch and turn off the power

- You cannot set personal data in the following conditions. Please re-register.
 - If no operation for more than 3 minutes when setting the data, the power will automatically switch off.
 - If the power goes off while you are setting.

Changing Personal Data

1

With the display unit inside the main unit, **Turn on the power**

omron 2012

After the display of current "year", "month and date", it will display "0.0 kg".

OMRON
REGGET

OMRON

Note: Do not step on the main unit or move the main unit before "0.0kg" is displayed. Do not put something on the main unit or touch display unit. Otherwise, it will display "Err (error)". (☑ page 27)

2 OMRON

Press the desired "Personal Profile Number" button.

The selected personal profile number lights.

3

Press DISPLAY button

"year of birth" flashes.

4

Changing personal data

Start from step 3 of "setting personal data". (page 11).

2 "0.0 kg" is displayed.

Personal data change is completed.

You can now take a measurement. (Fig. page 17)

5

Press power switch to turn off the power

Deleting Personal Data

1

OMRON

OMRON

OMRON

OMRON

OMRON

With the display unit inside the main unit,

Turn on the power

After the display of current "year", "month and date", it will display "0.0 kg".

Note: Do not step on the main unit or move the main unit before "0.0kg" is displayed, Do not put something on the main unit or touch display unit. Otherwise, it will display "Err (error)". (

2

Press "Personal Profile Number" button to select the desired number

The selected personal profile number lights.


3

Press DISPLAY button

"year of birth" flashes.


4

Deleting personal data

Press and hold the selected Personal Profile Number button more than 2 seconds

After "CIr" lights, return to set "date of birth".

Personal data is now deleted successfully.

5

Press power switch of the main unit and turn off power

Measuring Weight and Body Composition

(with Personal Profile Number Button / GUEST Button)

If you measure by using Personal Profile Number button, please preset personal data. (Fr page 11) If no personal data is available, use "GUEST" button to measure. (Fig. page 16)

Measure on flat and hard ground

Note: You cannot get correct measurement results on mat. carpet or other soft ground.


With the display unit inside the main unit, Turn on the power

After the displaying of current "year", "month and date", it will display "0.0 kg"

Note: • Do not step on the main unit or move the main unit before "0.0kg" is displayed. Do not put something on the main unit or touch display unit. Otherwise, it will display "Err (error)".

• When "year" and "month and date" of current date are displayed, it means the unit is ready to take a measurement. Do not move the main unit or put something on the main unit after power is turned on and before "0.0 kg" is displayed. When "Err" (error) is displayed, you will not get correct measurement results. (Frage 27)

Take out the display unit

Personal data have been set

Press "Personal Profile Number" button to select your own personal profile number.

The selected personal profile number lights.

■ When the right side screen appears...

There are no personal data under selected personal profile number. Set personal data at first. (Propage 11)

When measuring as a guest

What is "GUEST" mode?...

Users who does not set personal data in advance can measure weight, body fat and skeletal muscle by temporarily inputting "age, gender and height".

However, the unit will not display body index graph or memory value.

Press |GUEST| button

"G" lights and "age" flashes.

1 Set age

- Setting range is 10 ~ 99 years old.
- Press continuously to rapidly shift up or down.

After setting "age", "male / female" flashes.

Pollow the same steps to set "gender" and "height".


After setting "height", "age" and "gender", height will flash and 0.0kg will be displayed.

Measure weight and body composition (with Personal Profile Number button / GUEST button)

4

Taking a Measurement



Indicator flashes twice to show that weight has been determined.

Then measure body composition.

Keep staying on foot electrodes.


When "whole body" icon flashes,

Extend your arms straight at a 90° angle to your body


When measuring body fat and skeletal muscle, the same age comparison bar will display progress.


3 Display the measurement results


After measuring is completed, the same age comparison (whole body) will be displayed in accordance with weight, progress chart and different parts.

Step off the main unit.

Same age comparison at different parts (FF page 21 and 35)


Progress Chart (P page 24)

5 Check measurement results

Press $\binom{\text{DISPLAY}}{\text{SET}}$ button to switch display as shown in the following figure.


For viewing method of different parts, please press AREA CHECK button. "Viewing method of different positions display". (Propage 20)

Press button to view past measurement results.


Refer to "view previous measurement results (memory)" (

Viewing method of different positions display

The display changes to the next part with each press of the AREA CHECK button as follows:

Measure weight and body composition (with Personal Profile Number button / GUEST button)

6



Press power switch and turn off power

When viewing the previous measurement results after measurement is completed and power is turned off, refer to "Using the Memory Function" (Fig. 22).

Using the Memory Function

When measuring by using Personal Profile Number button, measurement results will be saved automatically. You can view measurement results of today, 1 day ago, 7 days ago, 30 days ago, 60 days ago and 90 days ago. When measuring with GUEST mode, measurement results cannot be recorded.

Note: • Personal profile number can record 97 days at most. After 97 days, the earliest measurement results will be deleted.

Measuring trace forgotten

Display 37 days ago

Measuring trace forgotten

Display 67 days ago

Measuring trace forgotten

Display 97 days ago

• If there is no past measurement results, " - - " will be displayed.


Measuring trace forgotten Measuring trace forgotten

Display 14 days ago

Display 6 days ago

Progress Chart

Indication bar of 5 stages will be used to display memory results of each index.


Final measurement results of "today" will be taken as normal value (0).

If you do not measure "today" but want to view past measurement results, it will be displayed with the latest memory results as beginning.

-If no measurement results exceed the age or measurement scope, indication bar will not display.

If today's measurement results exceed scope, progress Chart will not be displayed while memory results can be read.

You can confirm if reduction in weight is correct or not by viewing variation of other body conditions except weight.

As the example on the left, because of reduction in weight, the skeletal muscle increases. This is the correct method of reduction in weight but can not be detected easily.

Skeletal muscle (page 33)
What is correct reduction in weight? (page 38)

Location of Display Unit


■ The cord can rewind automatically.

Pull the display unit upwards to extend the cord attaching the display to the main unit. The cord is automatically rewound when you store the display in the main unit.

Pick up the display unit and rewind the cord into the unit.

If the cord does not rewind automatically...
 Untwist the cord, extend it once and try again.
 Do not extend the cord from the unit beyond the red marker line.

Push the display unit into the holder.

Measuring Weight only

Take a measurement directly without taking the display unit from the holder. It is unnecessary to set personal profile number.

Take a measurement on flat and hard ground

Note: You cannot get correct measurement results on mat, carpet or other soft ground.

1

With the display unit inside the main unit, **Turn on the power**

After the display of current "year", "month and date", it will display "0.0 kg".

Note: • Do not step on the main unit or move the main unit before "0.0kg" is displayed. Do not put something on the main unit or touch display unit. Otherwise measurement cannot

 If the current "year" and "month and date" displayed, please get ready for measuring.
 If you step onto the unit or put something before 0.0 kg appears on the display, an error message "Err" will appear. (Page 27) The measurement can not be performed correctly.

2

Step onto the unit

be performed correctly.


3

Check the measurement results

Indicator flashes twice to show that weight has been determined.

4

Press power switch, Turn off the power

Error Displays

■ Error Displays

Error Display	Cause	Correction
{r.	Your palms or soles are not in firm contact with the electrodes.	Press your palms or soles firmly to the electrodes, then measure. (🌃 page 17 and 18)
The posture for measurement is wrong or the palms or soles are not in firm contact with the electrodes. Measure without moving the hands or s (Firm page 17 and 18)		Measure without moving the hands or soles. (♣ page 17 and 18)
<pre>{rr}</pre>	The palms or soles are too dry.	Moisten the palms or soles with a wet towel, and then repeat measurement.
5113	Abnormal operation	Insert the batteries again and repeat measurement. Consult OMRON authorised dealer if problem persists.
	You stepped on the unit before the display indicated "0.0 kg".	Reconnect the power and step onto the main unit until the display indicated "0.0 kg".
C	You moved the main unit before the display indicated "0.0 kg".	Reconnect the power and do not move the main unit until the display indicates "0.0 kg".
Err	You moved your body during measuring.	Do not move your body during measuring.
	When measuring, body weight measuring range is exceeded.	Body weight of 135 kg or over is out of the measuring range of this unit. You cannot use this unit.

■ Troubleshooting

Problem	Cause	Correction	
"0.0 kg" is displayed even when you step onto the unit, or the weight value barely increases.	You stepped onto the unit before "0.0 kg" was displayed.	Wait until "0.0 kg" is displayed before stepping onto the unit.	
The value of body weight is displayed abnormally high or low.			
The value of body weight is displayed high	Display unit is taken out before "0.0 kg" is displayed.	Wait until "0.0 kg" is displayed before taking out display unit.	
	Measuring posture is incorrect	Take a measurement with correct posture (FFF page 17 and 18)	
	Is the unit placed on a carpet or other soft flooring? Is the floor surface is uneven?	Place the unit on a hard and level floor.	
Measurement results is abnormally high	The soles of your feet, or whole body are cold, causing poor blood circulation.	Warm yourself up, and then take a measurement once you feel your circulation is back to normal.	
or low Measurement results achieved each time s great different from each other	The electrodes are extremely cold.	Place the unit in a warm room for a few minutes, and then take a measurement when the electrodes no longer feel cold.	
	Palms and soles are too dry.	Slightly moisten your palms and soles with a damp towel and try again.	
" • • " is displayed for some of the results.	Your height, gender and birth is not correct This unit cannot measure body compositior if the height, gender and birth date settings	because of its measurement range, even	
	Batteries are not inserted.	Insert the batteries. (FF page 8)	
	The polarities of batteries are not aligned correctly.	Insert the batteries in correct alignment. (FF page 8)	
Nothing is displayed when the power is turned on.	The batteries are worn out.	Replace all four batteries with new ones. (page 8)	
	Cord connecting the main unit and the display unit is damaged.	Consult OMRON authorised dealer.	

■ Troubleshooting (continued from previous page)

Problem	Cause	Correction	
When weight measuring is	No personal data or GUEST identity is set (display unit does not show "personal profile number" or "GUEST")	Please set personal profile number or "GUEST" data before measuring again. (🏽 page 11 and 16)	
completed, weight and body composition cannot be measured.	"Personal Profile Number" button or button is pressed incorrectly (display unit does not show "personal profile number" or "GUEST")	Please press "Personal Profile Number" button or before take a measurement. (Page 16)	
Cord cannot be wound automatically	Cord is tangled.	Please pull out the cord and rewind. (I page 25)	
There is no operation and the power is turned off.	Please refer to "Automatic switch-off of the power" (FF page 28)		

X If measuring cannot be performed normally after the above has been done, please consult OMRON authorised dealer.

Automatic Switch-off of the Power

- If the monitor is not used within one minute of "0.0 kg" appearing on the display.
- If no information is entered for 3 minutes when entering personal data.
- If the monitor is not used for 3 minutes after "start measuring" or the measurement results are displayed.
- If the monitor is not used for 3 minutes after the target record is displayed.
- Personal profile number is not selected when measuring weight only. (about 30 seconds after body weight is displayed)
- When "Err" is displayed (about 10 seconds later)

X Although it is seldom seen, human physique difference might lead to incapability of measuring. In this case. Consult OMRON authorised dealer.

X If fault occurs or the unit is repaired, all personal data will be lost.

Specifications

Product	Body Composition Monitor HBF-375 (HBF-375-AP)			
	Display of Body Weight	0 to 100 kg with an increment of 100 g		
		100 to 135 kg with an increment of 200 g		
	ВМІ	2.5 to 90.0 with an increment of 0.1		
	Body Fat Percentage	5.0 to 50.0% with an increment of 0.1%		
	Body Fat Percentage Classification	- (Low) / 0 (Normal) / + (High) / ++ (Very High)		
	Segmental Subcutaneous Fat Per	centage (Whole Body, Trunk, Legs, Arms):		
		5.0 to 60.0% with an increment of 0.1%	Object age:	
	Segmental Skeletal Muscle Percentage			
	Whole Body	5.0 to 60.0% with an increment of 0.1%	years old	
Display*	Trunk	5.0 to 50.0% with an increment of 0.1%		
Diopiay	Legs	5.0 to 70.0% with an increment of 0.1%		
	Arms	5.0 to 60.0% with an increment of 0.1%		
	Resting Metabolism	385 to 3999 kcal with an increment of 1 kcal		
	Visceral Fat Level Visceral fat Classification	0.5 to 30 levels with an increment of 0.5 level 0 (Normal) / + (High) / ++ (Very High)	Object cas	
	Body Age	18 to 80 with an increment of one year	Object age: Over 17	
	Same Age Comparison	Subcutaneous Fat percentage: 5 levels	years old	
		Skeletal Muscle percentage: 5 levels	ľ	
	Progress Chart	5 levels		
	Memory	Today / 1 day ago / 7 days ago / 30 days ag	jo/	
	60 days ago / 90 days ago			
	Date	January 1, 2012 to December 31, 2030		
	Time 0:00 ~ 23:59 (24-hour display)			
	Setting Personal Data (the following	ng items can accommodate date of 4 people;)	
Set Items	Date of Birth	January 1, 1900 to December 31, 2030		
	Gender	Male / Female		
	Height 100.0 to 199.5 cm with an increment of 0.5 cm			
	Age (when using "GUEST" mode) 10 to 99 with an increment of one year * When using Guest mode, personal data should be set before measurement.			
	0.0 kg to 40.0 kg:	±400 g		
Weight Accuracy	40.1 kg to 135.0 kg:	±1%		
	4 AA batteries (R6)	2170		
Power Supply	(You may also use AA alkaline bat	teries (LR6).)		
Datta w. I ifa		ese batteries at room temperature of 23°C , n	neasured	
Battery Life	4 times a day)			
Operating				
Temperature /	+5℃ to +35℃ , 30% to 85% RH			
Humidity				
Storage Temperature /	-20°C to +60°C , 10% to 95% RH, 700 hPa - 1060 hPa			
Humidity / Air Pressure	Within well ventilated and dry room and no corrosive gas in surrounding air.			
Type of Protection Against Electric Shock	Type BF for applied part. Symbol: 1			
Weight	Approximately 2.6 kg (including batteries)			
External Dimensions	Approx. 329 (L) × 304 (W) × 54 mm (H)			
Accessory	Accessory 4 AA manganese batteries (R6), Instruction manual			

Subject to technical modification without prior notice.

^{*} For children under 10 years old: Only weight will be displayed.

For children over 9 and under 18 years old: All results except for Visceral Fat Level, Body Age and Same Age Comparison are displayed.

For people 81 years old or older: measured body composition can only be used for reference.

For 18 - 80 years old, all results are displayed.

Information on Body Composition

BMI & IDEAL WEIGHT

Everybody is concerned about their weight. But what is the basis for determining your ideal weight?

The key is BMI.

BMI is short for Body Mass Index. This index uses the following simple formula to indicate the ratio between weight and height of a person.

BMI = weight (kg) / height (m) / height (m)

The OMRON Body composition monitor uses the Height information stored in your personal profile number to calculate your BMI.

The Ideal BMI is 22

Maintaining an ideal weight can help prevent obesity or weight loss and other diseases, and lead a longer life. The ideal weight for a BMI of 22 is calculated as follows.

Ideal Weight (kg) = $22 \times \text{height (m)} \times \text{height (m)}$

Example: For a person 1.55 m tall...

 $22 \times 1.55 \times 1.55 = 52.855 \rightarrow \text{Approx } 52.9 \text{ kg}$

However, this method of ideal weight calculation may not be applicable for professional athletes and body builders, who have higher muscles ratio in their bodies.

Interpreting the BMI Result

BMI Classifications (by the Wh		
Less than 18.5	Underweight	
18.5 or more and less than 25	Normal	
25 or more and less than 30	Overweight	
30 or more	Obese	

The above-mentioned indices refer to the values for obesity judgment proposed by the WHO, the World Health Organization.

INFORMATION ON BODY COMPOSITION

BODY FAT

Body fat serves a vital role in storing energy and protecting internal organs. We carry two types of fat in our bodies: 1) essential fat which is stored in small amounts to protect the body and 2) adipose tissue (stored fat) provides cushion and insulation to internal organs, covers the nerves, moves vitamins throughout the body and is the largest reserve of stored energy available for activity. While too much body fat may be unhealthy, having too little fat can be just as unhealthy. Also, the distribution of body fat in men and women is different, so the basis for classifying the body fat percentage is different between the genders.


What is Body Fat Percentage?

Body fat percentage refers to the amount of body fat mass in regards to the total body weight expressed as a percentage.

Body fat percentage (%) = {Body fat mass (kg) / Body weight (kg)} \times 100 Our device uses BI method to estimate your body fat percentage.

Depending on where fat is distributed in the body, it is classified as visceral fat or subcutaneous fat

Interpreting the Body Fat Percentage Result


Based on the obesity values proposed by Lohman (1986) and Nagamine (1972).

VISCERAL FAT

Visceral fat is found in the abdomen and surrounding vital organs. It is different from fat found directly underneath the skin, which is referred to as subcutaneous fat. Visceral fat can go largely unnoticed because it's not visible to the naked eye. One way visceral fat can be seen is through Magnetic Resonance Imaging (MRI).

MRI Cross Section at the Navel Level

Too much visceral fat is thought to be closely linked to increased levels of fat in the bloodstream, which may lead to conditions such as high cholesterol, heart disease and type 2 diabetes. In order to prevent or improve these conditions, it is important to try to reduce the amount of visceral fat levels to an acceptable level.

Interpreting the Visceral Fat Level Result

Visceral Fat Level Level Classification		
0.5 - 9.5	0 (Normal)	
10.0 - 14.5	+ (High)	
15.0 - 30.0	++ (Very High)	

^{*} The data referenced by OMRON Healthcare.

NOTE: Visceral fat levels are relative and not absolute values.

What is Subcutaneous fat?

Subcutaneous fat = fat below the skin

Subcutaneous fat not only accumulates around the stomach but also around the upper arms, hips and thighs, and can cause a distortion of the body's proportions. Although not directly linked to increased risk of disease, it is thought to increase pressure on the heart and other complications.

RESTING METABOLISM

Regardless of your activity level, a minimum level of energy is required to sustain the body's everyday functions. Resting metabolism, the amount of calories needed to supply the body with the minimum level of energy, differs between individuals depending on variables such as age, weight, body composition, and energy expenditure.


60-70% of daily energy use is for resting metabolism

The total amount of energy used by the body in a typical day is as follows:

Resting metabolism	Energy required to maintain vital functions.
Daily activity metabolism	Energy used for daily activities such as commuting to work, household chores, hobbies etc.
Diet-induced thermogenesis	Energy emitted after eating a meal.

The ratio of these is 60%-70% for resting metabolism, 20%-30% for daily activity, and 10% for diet induced thermogenesis. This means that resting metabolism accounts for most of our daily energy consumption.

If our daily caloric intake exceeds the amount of energy required for these activities, the additional calories can be stored as fat.

SKELETAL MUSCLE

What is skeletal muscle ... Muscle is divided into two types, muscle in internal organs, such as the heart, and muscle attached to bones that is used to move the body. Skeletal muscle can be increased through exercise and other activity.

What is skeletal muscle ratio · · · Refers to the amount of skeletal muscle mass in regards to the total body weight expressed as a percentage.

Skeletal muscle ratio (%) = {Skeletal muscle mass (kg) / Body weight (kg)} \times 100

This product uses BI (Bioelectrical Impedance) to calculate your skeletal muscle mass.

Keep skeletal muscle and strengthen it to make corporeity not easily become fat.

This product feature is the measurement of skeletal muscle ratio. If you strengthen skeletal muscle and improve basal metabolism, you will not get fat easily.

Interpreting the Skeletal Muscle Percentage Result (for adults)

Gender	– (Low)	0 (Normal)	+ (High)	++ (Very High)
Female	5.0-25.8	25.9-27.9	28.0-29.0	29.1-60.0
Male	5.0-32.8	32.9-35.7	35.8-37.3	37.4-60.0

According to OMRON HEALTHCARE figures

BODY AGE

Body age is based on your resting metabolism. Body age is calculated by using your weight and body fat percentage to produce a guide to judge whether your body age is above or below the average for your actual age. To know your body age is useful to improve your health condition.


Body age varies according to Body composition and resting metabolism, even if your height and weight is the same

The figures below shows some examples of body age.

The example on the left has a body age of 30, or average body composition for the actual age.

The example in the middle has a lower body fat percentage than average and higher resting metabolism, which results in a body age of 25 that is younger than the actual age.

The example on the right has a higher body fat percentage than average and lower resting metabolism, which results in a body age of 35 that is older than the actual age.


※ For children under 18 years old: Visceral fat level, body age and same age comparison will not be displayed.

Same Age Comparison

Understand body condition through same age comparison

Subcutaneous fat and skeletal muscle at different parts are displayed. Its level is displayed through comparison with same age and same BMI.

The average is investigation result of our company. Please do not judge health condition in accordance with it.

^{*} The monitor cannot display same age comparison for users under 18 years old.

Feature and Measuring Principle of Body Composition Monitor

Feature of body composition monitor

Is the method of whole body measuring, which is not subject to daily variation.

1 Water content in our body tends to accumulate in our leg during day. That's why we usually feel dropsy at the feet at dusk or in the evening.

In the morning and dusk, body resistance distribution will change, as water proportion in upper body and lower body changes.

KaradaScan measures whole body resistance by using electrode of both hands and both feet. Therefore, it is less subject to variation of water content distribution in human body.

Basic data used to establish body composition formula, Collect by using measuring method and measuring device most suitable for various body compositions.

Body fat: weighing method in water (Submerge your head into water pool and exhaust all air in your lung before measuring body weight. The method to calculate body fat in accordance with body density is based on Archimedes principle.)

Resting metabolism: breath analyzing device (used to analyze air composition in breath). It can be used to measure oxygen consumption of breath at quiet time to calculate resting metabolism.)

Subcutaneous fat and skeletal muscle: DXA (Dual Energy X-Ray Absorptiometry) method (measuring body composition at different parts by quantifying bone quantity and body composition of 2 kinds of X rays with different wave length without the need of resolving body to calculate weight.)

MRI (Magnetic Resonance Imaging System) (large number of hydrogen atomic nucleus contained in human body resonate with magnetic field and electric wave. Image hydrogen atomic nucleus and nuclear magnetic resonance to calculate human body composition.)

Measuring principle of body composition monitor

Calculate body fat in accordance with resistance value.

KaradaScan measures body resistance by using weak currency flowing through both hands and both feet (BI method*).

* BI method (Bioelectrical Impedance/ biological resistance method).

Tissues with more water content in human body tend to conduct electricity (such as muscle and vein) easily. Fat tissue almost conducts no electricity. The body feature is used to calculate rate of fat tissue and non-fat tissue.

Currency flowing through human body is very weak (50 kHz, $500\mu A$), which is not stimulant and is very safe to human body.

In order to obtain body fat and other data from resistance between both hands and both feet,

The five items, i.e., (1) resistance value, (2) height, (3) weight, (4) age and (5) gender are required, which are obtained in accordance with basic human data collected by the company independently.

Things to Take Note

Proper measuring place... flat and hard ground

You cannot get correct measurement results on mat, carpet or other soft ground.

Correct posture... arm and body are kept at 90°

In order to obtain resistance value necessary to calculate body fat, you should basically keep your arm and body at 90° and straighten your arms. In order to ensure stability of measurement results, please take a correct posture when performing measuring.


Suitable period for measuring... After getting up and 2 hours after taking food

Water content of the part where electric flows through and which is measured by the product will affect body resistance.

Please refer to the following figure and perform measuring while water content in body is steady.

Recommend Measurement Time period

In case of any doubt about measurement results... (I T in case of doubt in instruction manual, page 27)

Do not perform measuring in the following cases.

- After strenuous exercise
- After sauna or shower
- After drinking much alcohol
- After taking in much water or food (1~2 hours)


Note: The following people may not accurately measure their body fat and other index, as water content in their body and body composition are greatly deviated from average value.

- Children in growth stage
- People of high age
- Patients with fever due to cold
- Pregnant women
- Patients having osteoporosis
- Patients having dropsy
- Patients receiving hemodialysis
- Professional sportsmen or body builders (or similar people)
- * It is recommended that above stated people confirm "variation tendency" of their body fat and other index before using the product.

Q&A of Body Composition Monitor

■ What is correct weight reduction?

If you lose weight by going on a diet instead of taking exercise and neglect nutrition balance, even if your weight is reduced, resting metabolism will decrease as muscle (skeletal muscle) decreases and you are more likely to become fatter.

In order to avoid reversing,

You should build up non-fat physique by increasing skeletal muscle and improving resting metabolism.

Extreme weight losing is most likely to reverse your body weight. When body weight is reversed, visceral fat is more likely to accumulate in human body than subcutaneous fat. Visceral fat is the major cause of life habit related disease. It is body weight reversion that leads to visceral fat obesity.


Pay attention to visceral fat level (**№** page 32)

■ If you want to be a person not tending to be fat...

You should lay emphasis on resting metabolism (Page 33)

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OMRON HEALTHCARE Co., Ltd.
OMRON HEALTHCARE SINGAPORE PTE LTD.