
Digital Weight Scale
Model HN-286

- Thank you for purchasing Omron Digital Weight Scale HN-286. Be sure to read the instruction manual before use and practice the instructions.
- For safe use of this product, please read the "Safety Precautions" on the Instruction Manual firstly.
- Please take good care of the Instruction Manual for future reference.

Notes on Safety

- The warning sign and legend shown here are to ensure proper and safe use of products, and to prevent harm or damage to the user or others.
- Warning sign and its meaning shown as follows.

Warning sign	Meaning
 Warning	Casualties may occur during wrong use.
 Caution	Personal injury or article damage may occur during wrong use.

Legend	
 Disassembling prohibited	⊘ operation prohibited or not allowed. Prohibited matters are indicated inside or around ⊘ in forms of text or picture. The left picture means "Disassembly prohibited".
 Mandatory:	● mandatory (rules must be followed). Mandatory matters are indicated inside or around ● in forms of text or picture. The left picture means "generally mandatory issues".

Warning	
● Do not step on the edge of the scale. • You may fall and get injured, and accuracy cannot be guaranteed.	 Prohibited
● Do not put on wet floors, including tiles or other slippery surfaces. • You may slip and get injured.	 Prohibited
● In case you are in a diet or exercise program for weight loss, do not judge by yourself. Please follow the guidance of a doctor or specialist. • Self-judgement may be detrimental to your health.	 Mandatory
● Please use indoor.	 Mandatory

Warning	
● Do not step on the scale with wet feet. • You may slip and get injured.	 Mandatory
● Do not jump to the scale or bounce on it. • You may slip and get injured. • The scale may be damaged by the impact of weight beyond specified scope.	 Prohibited
● Do not put your fingers into the gap or the inside of the scale. • You may get injured.	 Prohibited
● Do not disassemble the scale. • The scale may be damaged during disassembling.	 Disassembling prohibited
● If battery fluid gets into eyes or touches skin, clothing, rinse with plenty of water immediately. • It may lead to blindness or injury. Please seek professional medical treatment.	 Mandatory
● Do not throw the battery into fire. • It may rupture and cause injuries.	 Prohibited

Caution	
● Please place it on a solid flat surface before use. • Or else you may slip and get injured, and cause measurement errors.	 Mandatory
● Do not place the scale upright against the wall or heavy objects on it, as measurement error or unit damage may occur. ● Please place it beyond the reach of infants and young children to avoid danger. ● When used by customers with tinea pedis or other skin diseases, it shall be disinfected to avoid cross infection. ● The scale is for home use only. ● Do not use mobile phone near the scale, which may cause inaccurate measurement. ● Do not use single hand to hold the unit because its glass surface is smooth.	 Prohibited
● Do not put in damp and moist environment, under direct sunlight, or near air conditioning or fires. • Malfunction may ensue.	 Prohibited
● Do not wash the scale with water. • It is likely to cause faults.	 Prohibited

● The weight scale is for home use. • Do not use for purposes other than measuring body weight.	 Prohibited
● People with disabilities, or who are physically frail, should always be assisted by another person when using this unit. Use a handrail or so when stepping on the unit.	 Mandatory
● It may not measure correctly on soft ground such as on a carpet. Please put it on a solid and flat surface.	 Mandatory

Tips

● Be sure to use the specified type of batteries.	
● When installing the battery, do not mistake the directions of ⊕ ⊖ pole.	
● During out-of-service period (about 3 months or above) of the product, remove the battery.	
● When discarding the battery, please follow relevant provisions on environmental protection.	

Product Features

- (1) High precision: the smallest display unit is 0.1Kg.
- (2) Large-size Display: 38mm x 77mm
- (3) High standard toughened glass is used on the surface of the weight scale.
- (4) Automatic startup feature (for immediate measurement), and automatic power off features (power is cut off 16 seconds after measurement).

About the Accuracy

- | | | |
|--|---------------|--|
| (1) Rigorous test will be exercised before delivery. The accuracy of a newly produced scale ensured by the manufacturer is as shown in the right table. | The weight | The accuracy of a newly produced scale ensured by the manufacturer |
| | 5kg-60kg | ±400g |
| | 60.1kg-100kg | ±600g |
| (2) The scale is manufactured for measuring body weight.
Note: According to the Measurement Law, the product can not be used for measuring tradable goods or be used as other weighting apparatus for test purpose. | 100.1kg-180kg | ±1% |
- (3) The measurement range of the scale is 5kg-180kg.

Correct Method of use

Battery placement

Please put in one 3V lithium battery before use.
The scale depends on the attached battery. Please pull out the battery insulation film before use and change the battery when it low power.

- 1 Open the battery cover.
- 2 Put one 3V lithium battery into the battery box (positive pole upward) and then tightly close the battery box.

When you install the battery for the first time

1	Full display after battery placement	
2	Display after 1-2 seconds Make initial setup (zero adjustment before measurement)	
3	Place the scale horizontally on a flat ground. The power turns off after the scale displays the value automatically (Note: the value 1.5 is for instance) Note: If it is not placed horizontally, incorrect measurement may occur.	


Measurement

The scale will show the value when the user stands stably on the center of platform.
When (kg) is displayed, the measurement ends. The scale will be automatically turned off after 16 seconds. If you stand in an unstable manner, the display value will flash and automatically shut down 10 seconds later.

Error Display

Error symbol	Cause	Solutions
	The battery is running out	Change the battery
	The weight measured is beyond the measurement scope	
	Measurement error	Repeat Measurement

About the battery

- (1) If the display shows  it means the battery is running out. Please replace the battery.
- (2) The battery is one 3V lithium battery.
- (3)  sign may not appear when old battery is used.
- (4) Battery life is approximately 1 year (measured 4 times a day).
- (5) Remove the battery when not using the unit for a long time. In addition the unit with installed battery is not used for long time, it may cause battery leakage and product failure.

Maintenance and Storage

- Do not use gasoline, thinner or hot water to clean it. Use a soft cloth with mild detergent to carefully wipe the stains on the scale and then wipe it with another soft cloth.
- If an item is placed on the weight scale for a long time, it may result in inaccurate measurements or damage the scale.
- Please don't place it upright because that result inaccurate measurements.
- In order to measure more accurately, this product has the function of auto zero adjustment. If the measurement result is inaccurate due to placing heavy objects on it or placing it upright against the wall, this scale can carry through zero adjustment automatically after every measurement, please measure again.

Troubleshooting

Problem	What to check	Corrections
Fail to turn on	Is the scale placed on uneven ground?	Please measure on hard level ground.
	Is the battery running out?	Change the battery.
	Is battery polarity ⊕ ⊖ placed correctly?	Please insert ⊕ ⊖ battery according to the right direction.
Weight showed is unusually light or heavy	Is the scale placed on uneven ground?	Please measure on hard level ground.

Specifications


Product Name	Digital Weight Scale
Model	HN-286
Weight Range	5kg - 180kg
Minimum Displaying Unit	0.1kg
Power	1 battery 3V lithium battery (CR2032 type)
Operating Temperature and Humidity	+5°C-+35°C/ 30%RH-85%RH
Storage Temperature and Humidity	-20°C-+60°C/ 10%RH-95%RH
Atmospheric Pressure for Operation	860hPa-1060hPa
Atmospheric Pressure for Transportation and Storage	700hPa-1060hPa
Battery Life	Approximately 1 year (measured 4 times a day)
Dimensions	Approximately (L)300mm x (W)300mm x (H)18mm
Weight	Approximately 1.5Kg (including battery)
Accessories	Instruction manual, a trial battery

※ Specifications are subject to change without prior notice.

	OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, Kyoto 617-0002 JAPAN
Asia Pacific HQ	OMRON HEALTHCARE SINGAPORE PTE LTD. 438A Alexandra Road, #05-05/08 Alexandra Technopark, Singapore 119967 www.omron-healthcare.com.sg

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